ENGLISH

INSTRUCTIONS FOR USE
MIVI Q™ Aspiration Catheter

Contents
1X Catheter with pin vise.

Device Description

The MIVI Q Aspiration Catheter is a single-lumen, variable stiffness catheter with radio-paque markers on the distal and proximal end of the catheter portion for angiographic visualization. The catheter shaft has a hydrophilic coating to reduce friction during use. The proximal portion of the catheter is a stainless-steel control wire. The Q Aspiration Catheter may be introduced via an 8F guide catheter / 6F guide sheath and over a guidewire / microcatheter into the arterial vasculature until the desired vessel is reached. A vacuum source is attached to the guide catheter/sheath and aspiration is performed from the Q Aspiration Catheter through the guide/sheath to aspirate thrombus.

Device Dimensions

Device Model	Length (cm)	Proximal Catheter Section ID	Proximal Catheter Section OD	Distal Catheter Section ID	Distal Catheter Section OD	Cath-eter Sec-tion Length (cm)
Q6-69145-E	145	1.75 mm (0.069")	2.24 mm (0.088") 6.7F	1.75 mm (0.069")	2.13 mm (0.084") 6.4F	25
Q5-57145-E	145	1.45 mm (0.057")	2.24 mm (0.088") 6.7F	1.45 mm (0.057")	1.83 mm (0.072") 5.5F	25
Q4-43150-E	150	1.45 mm (0.057")	2.24 mm (0.088") 6.7F	1.09 mm (0.043")	1.40 mm (0.055") 4.2F	30
Q3-36163-E	163	1.45 mm (0.057")	2.24 mm (0.088") 6.7F	0.91 mm (0.036")	1.22 mm (0.048") 3.7F	43

Device Compatibility

Device Description	Maximum Guidewire Size	Maximum Microcatheter Size	8F Guide Catheter/ 6F Sheath ID Range
Q6-69145-E	0.035"	5.0F	0.088" – 0.090" (2.2 – 2.3 mm)
Q5-57145-E	0.035"	4.0F	0.088" – 0.090" (2.2 – 2.3 mm)
Q4-43150-E	0.035"	3.0F	0.088" – 0.090" (2.2 – 2.3 mm)
Q3-36163-E	0.018"	2.4F	0.088" – 0.090" (2.2 – 2.3 mm)

Indications for Use

The Q Aspiration Catheter is indicated for the removal of fresh, soft emboli and thrombi in the peripheral and neurovascular systems. It may also be used as a diagnostic angiographic catheter.

Contraindications

• Vessels smaller than the catheter’s distal outer diameter.

WARNINGS

• The Q Aspiration Catheter is intended to be used by physicians skilled in the field(s) of Interventional Radiology, Interventional Neuroradiology, Neurosurgery and/or Neurology who are skilled in endovascular interventional procedures for acute ischemic stroke or acute limb ischemia.
• Do not use the device in patients who are allergic to polyether block amides (PEBAX), polytetrafluoroethylene (PTFE) or polyether based thermoplastic polyure-thanes (TPU).

• Never advance or retract the catheter against resistance without careful assess-ment of cause using fluoroscopy. If cause cannot be determined, withdraw the catheter.
• Do not torque the control wire of the Q Aspiration Catheter. Torquing, movement against resistance or forced insertion of the catheter may result in device damage, vessel damage and/or tip breakage which may result in injury to the patient.
• If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove or replace catheter.

Precautions

• Store in a dark, dry place.
• Use prior to the "Use By" date.
• Do not expose device to solvents.
• The device is intended for single use only. Do not resterilize or reuse. Resteriliza-tion and/or reuse may result in impaired structural integrity and/or function.
• Do not use open or damaged packages.
• Upon removal from package, inspect device to ensure it is not damaged. Do not use kinked or damaged devices.
• Confirm the compatibility of the Q Aspiration Catheter with any interventional devices to be used in conjunction with the catheter (e.g. guide catheter, micro-catheter, guidewire).
• Exercise care when handling the catheter before and during a procedure to reduce the possibility of damage.
• Once hydrated do not allow the catheter to dry out.
• Upon insertion and removal of the catheter, ensure the hemostasis valve is large enough to accommodate smooth movement.
• Manipulate the catheter under fluoroscopic visualization. Do not attempt to move the catheter without observing the resultant tip response.
• Avoid angular strain on the control wire during preparation and insertion.
• Ensure that the proximal radiopaque marker does not extend past the distal marker band of the guide catheter.
• When performing aspiration, ensure that the stopcock on the rotating hemostasis Y valve is open for only the minimum time needed to remove thrombus. Excess aspiration or failure to close the stopcock when aspiration is complete is not recommended.
• Excessive force on the Q Aspiration Catheter is not recommended. If repositioning of the catheter is necessary, repositioning should be performed using appropriate techniques.
• Do not exceed 300 kPa during contrast injection.
• Do not use automated high-pressure contrast injection equipment with the Q Aspiration Catheter because it may damage the device.
• If an intraluminal device becomes lodged in catheter, or if the catheter becomes severely kinked, withdraw the entire system (intraluminal device, catheter and guiding catheter).

Potential Complications

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with possible complications.

Potential complications include, but are not limited to, the following:

• access site complications (hemorrhage, hematoma, fistula, etc.)
• air embolism
• allergic reaction (contrast, materials, device)
• death
• device malfunction
• emboli/distal embolization
• inability to remove all thrombus
• infection
• intracranial hemorrhage
• ischemia
• kidney damage
• neurological deficits or new stroke
• vasospasm
• vessel trauma (dissection or perforation)

Devices Required for Use with the Q Aspiration Catheter

• 8F arterial guiding catheter or a 6F long femoral sheath
• Rotating hemostasis valve with Y adapter and stopcock
• Compatible guidewire (recommended)
• Compatible microcatheter (recommended)
• Appropriate vacuum source for aspiration

Preparation and Directions for Use

1. Open the shelf box and inspect pouch for damage.
2. Open pouch using standard sterile handling procedure.
3. Remove catheter from package and inspect for kinks or damage. Do not use a kinked/damaged device.
4. Hydrate with saline before use. Once hydrated, do not allow the Q Aspiration Catheter to dry out.
5. Introduce Q Aspiration Catheter with microcatheter into vasculature through a guide catheter/introducer sheath and over an indwelling guidewire using percutaneous access technique of choice.
6. Under fluoroscopic guidance, advance the Q Aspiration Catheter over the microcathe-ter/guidewire until desired position is attained.
7. Carefully remove guidewire and microcatheter.
8. Advance the tip of the Q Aspiration Catheter to the thrombus.
9. Attach an appropriate vacuum source to stopcock of the guide catheter/ sheath.
10. Activate vacuum source to begin aspiration.
11. Once aspirated fluid is visible, position the catheter as needed to capture target thrombus.
12. When aspiration is complete, stop the vacuum source.
13. If there is clot remaining, flush catheter and repeat steps above for a maximum of three passes.
14. Carefully remove and discard all devices at the end of the procedure.

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